



Comfiwave[®]

Compression to drift off in

What is *Comfiwave*®?

Comfiwave is a patented **compression garment**, designed to complement the treatment a patient receives during waking hours.

Comfiwave gives a therapeutic pressure of **15-21mmHg**.

Easy donning and doffing gives Comfiwave an innovative versatility, providing both night-time compression and a comfortable solution when resting in the daytime.



Below Knee
Page 10 ►



Thigh High
Page 12 ►



Toecap
Page 14 ►



Shorts
Page 16 ►



Sleeve
Page 18 ►



Sleeve + Mitten
Page 20 ►



Glove
Page 22 ►



Glove to Elbow
Page 24 ►



Glove to Axilla
Page 26 ►




Abdominal
Page 28 ►



Breast Band
Page 30 ►



Neck
Page 32 ►



Kids
Page 34 ►



For night time

Oedema related swelling doesn't stop overnight but traditional compression garments are not best suited to being worn in bed. Continue compression therapy overnight in supreme comfort with Comfiwave's soft and breathable fabric providing gentle compression while you sleep.



For any time

Comfiwave garments use proprietary 'low-resistance' yarns which allow them to be easily donned and doffed. Perfect for relaxing during the day.

The breathability of the fabric means that Comfiwave is suitable for use all year round, being cooler to wear than conventional hosiery.

Whole-garment knitting technology

Easy colour coded sizing system

e.g. ■■■ = Large, Long

Comfiwave's wave patterning mimics the pressure differential massage effect that is ordinarily found in cut-foam, chip-bag style garments.

Comfiwave lower limb garments utilise a stiffened and shaped ankle to prevent 'cutting in' of the anterior ankle

The lateral toe section has also been knitted to stretch more, which can prevent pressure sores on the outside of the toes

Heel section helps prevent slippage



Comfiwave® features



Easy to don and doff

Patented, high-elastic knit makes Comfiwave easy to put on and take off



Made from a soft, undyed cotton

Soft cotton against the skin and seamless design ensures maximum comfort



Breathable fabric

Comfiwave's wave patterning gives it high air permeability for coolness and stretch

Comfiwave[®] indications



Provides safe, night-time compression

Complements treatment a patient receives during waking hours



Can be used in palliative care



Can be worn under easywrap or bandaging



I love this garment. I have been wearing it overnight for the past week and it has made a massive difference!

- Lymphoedema patient

Really comfortable to wear at night. Unlike other products doesn't make my leg hot when temperature is high. Easy to put on and remove. Excellent product.

- Lymphoedema patient

“Patients reported that Comfiwave was comfortable and easy to use, which itself would improve concordance with treatment and self-management of oedema.”

“...this type of product can be used for stand-alone treatment or as an adjunct to other treatment modalities and for 24 hours a day if needed.”

Download
clinical case
studies and article

www.comfiwave.info


*British Journal of
Community Nursing*



Comfiwave® Below Knee

	SHORT		REGULAR		LONG
S	CW1-AD2S 422-1412	■	CW1-AD2R 422-1446	■	CW1-AD2L 422-1438
M	CW1-AD3S 422-1362	■	CW1-AD3R 422-1420	■	CW1-AD3L 422-1479
L	CW1-AD4S 422-1370	■	CW1-AD4R 422-1396	■	CW1-AD4L 422-1453
XL	CW1-AD5S 422-1388	■	CW1-AD5R 422-1404	■	CW1-AD5L 422-1461

Drug Tariff codes shown in red



Circumferences cm

S

M

L

XL

c **Widest Calf**

25 - 32

30 - 38

36 - 46

44 - 56

b **Ankle**

17 - 22

20 - 26

22 - 30

27 - 35

Lengths cm

SHORT

REGULAR

LONG

A-d **Floor to
Below Knee**

32 - 37

37 - 42

42 - 48



b
Ankle



c
Widest Calf




A-d
Floor to below knee



Comfiwave® Thigh High

	SHORT		REGULAR		LONG	
S	CW1-AG2S 422-1487	■	CW1-AG2R 422-1529	■	CW1-AG2L 422-1586	■
M	CW1-AG3S 422-1560	■	CW1-AG3R 422-1537	■	CW1-AG3L 422-1495	■
L	CW1-AG4S 422-1552	■	CW1-AG4R 422-1511	■	CW1-AG4L 422-1545	■
XL	CW1-AG5S 422-1503	■	CW1-AG5R 422-1578	■	CW1-AG5L 422-1594	■
XXL	CW1-AG6S 429-5507	■	CW1-AG6R 429-5515	■	CW1-AG6L 429-5499	■

Drug Tariff codes shown in red



Circumferences cm

S

M

L

XL

XXL

g Thigh Top

40 - 48

46 - 56

54 - 65

62 - 75

72 - 92

c Widest Calf

25 - 32

30 - 38

36 - 46

44 - 56

52 - 64

b Ankle

17 - 22

20 - 26

22 - 30

27 - 35

27 - 35

Lengths cm

SHORT

REGULAR

LONG

A-g Floor to Thigh Top

64 - 70

70 - 77

77 - 84




b
Ankle



c
Widest Calf



g
Thigh Top



A-g
Floor to thigh top



Comfiwave® Toecap

REGULAR

M

CW1-ZA3



L

CW1-ZA4





Circumferences cm

M

L

a Ball of Foot

20 - 27

26 - 34




a
Ball of foot



Comfiwave® Shorts

REGULAR

S	CW1-GT2	
M	CW1-GT3	
L	CW1-GT4	
XL	CW1-GT5	
XXL	CW1-GT6	



Circumferences cm

S

M

L

XL

XXL

t

WAIST

50 - 68

66 - 86

84 - 104

102 - 124

120 - 146

g

THIGH TOP

40 - 48

46 - 56

54 - 65

62 - 75

72 - 92




g

Thigh Top

Comfiwave® Sleeve

	SHORT		REGULAR		LONG	
S	CW1-CG2S 422-1198	■	CW1-CG2R 422-1172	■	CW1-CG2L 441-8018	■
M	CW1-CG3S 422-1180	■	CW1-CG3R 422-1149	■	CW1-CG3L 441-8026	■
L	CW1-CG4S 422-1164	■	CW1-CG4R 422-1156	■	CW1-CG4L 441-8034	■

Drug Tariff codes shown in red



Circumferences cm		S	M	L
g	Axilla	25 - 30	30 - 40	40 - 50
d	Mid-point of forearm	23 - 26	26 - 30	30 - 38
c	Wrist	13 - 18	18 - 22	22 - 28

Lengths cm		SHORT	REGULAR	LONG
c-g	Wrist to Axilla	36 - 42	40 - 46	44 - 50



c
Wrist




d
Mid-point of
forearm



g
Axilla




c-g
Wrist to axilla



Comfiwave® Sleeve + Mitten

	SHORT		REGULAR	
S	CW1-AH2S 422-1214	■	CW1-AH2R 422-1230	■
M	CW1-AH3S 422-1206	■	CW1-AH3R 422-1248	■
L	CW1-AH4S 422-1222	■	CW1-AH4R 422-1255	■

Drug Tariff codes shown in red



Circumferences cm

S

M

L

g Axilla

25 - 30

30 - 40

40 - 50

d Mid-point of forearm

23 - 26

26 - 30

30 - 38

c Wrist

13 - 18

18 - 22

22 - 28

a Pal at base of fingers

16 - 22

18 - 24

20 - 26

Lengths cm

SHORT

REGULAR

c-g Wrist to axilla

36 - 42

40 - 46



a

Palm at base
of fingers



c

Wrist



d

Mid-point of
forearm



g

Axilla



c-g

Wrist to axilla



Comfiwave® Glove

REGULAR

S

CW1-ZC2



M

CW1-ZC3



L

CW1-ZC4





Circumferences cm		S	M	L
b	Base of thumb	17 - 22	20 - 25	23 - 28
c	Wrist	14 - 20	16 - 24	22 - 28



b
Base of thumb
webbing



c
Wrist



Comfiwave® Glove to Elbow

REGULAR

S

CW1-ZE2



M


CW1-ZE3



L

CW1-ZE4





Circumferences cm

b Base of thumb webbing

S

17 - 22

M

20 - 25

L

23 - 28

c Wrist

14 - 20

16 - 24

22 - 28

d Mid-point of forearm

23 - 28

26 - 32

30 - 38

Lengths cm

REGULAR

c-e Wrist to elbow

20 - 24



b
Base of thumb
webbing




c
Wrist



d
Mid-point of
forearm



c-e
Wrist to elbow



Comfiwave®

Glove to Axilla

SHORT

REGULAR

S

CW1-ZH2S

CW1-ZH2R

M


CW1-ZH3S

CW1-ZH3R

L

CW1-ZH4S

CW1-ZH4R



Circumferences cm		S	M	L
g	Axilla	25 - 30	30 - 40	40 - 50
d	Mid-point of forearm	23 - 26	26 - 30	30 - 38
c	Wrist	13 - 18	18 - 22	22 - 28
a	Palm at base of fingers	16 - 22	18 - 24	20 - 26
Lengths cm		SHORT	REGULAR	
c-g	Wrist to Axilla	36 - 42	40 - 46	



a
Palm at base
of fingers



c
Wrist



d
Mid-point of
forearm



g
Axilla



c-g
Wrist to axilla




Comfiwave® Abdominal Band

M	CW1-HT3 436-1556	■ ■ ■
L	CW1-HT4 436-1564	■ ■ ■
XL	CW1-HT5 436-1549	■ ■ ■



Comfiwave® Thigh + Hip

	SHORT		REGULAR	
S	CW1-ET2S	■ ■	CW1-ET2R	■ ■
M	CW1-ET3S	■ ■	CW1-ET3R	■ ■
L	CW1-ET4S	■ ■	CW1-ET4R	■ ■
XL	CW1-ET5S	■ ■	CW1-ET5R	■ ■



Circumferences cm

t Waist

M

60 - 75

L

70 - 95

XL


90 - 120

j Hips

80 - 100

90 - 115

110 - 140



Circumferences cm

g Thigh top

40 - 48

46 - 56

54 - 65

62 - 75

e Knee

27 - 35

33 - 39

41 - 50

48 - 60

Lengths cm

SHORT

REGULAR

e-g Knee to thigh top

20 - 26


25 - 31





Comfiwave® Breast Band

S	CW1-TN2 436-1655	■ ■ ■ ■
M	CW1-TN3 436-1648	■ ■ ■ ■
L	CW1-TN4 436-1630	■ ■ ■ ■
XL	CW1-TN5 436-1770	■ ■ ■ ■



Circumferences cm

Above breast

S

70 - 88

M

80 - 102

L

90 - 116

XL

100 - 130

Over breast

70 - 92

80 - 108

90 - 126

100 - 144

Below breast

70 - 80

80 - 92

90 - 104


100 - 116



Above Breast



Over Breast




Under Breast



Comfiwave® Neck

	SHORT		REGULAR	
S	CW1-HN2S	■	CW1-HN2R	■
M	CW1-HN3S	■	CW1-HN3R	■
L	CW1-HN4S	■	CW1-HN4R	■



BASE OF NECK

S

32 - 38

M

37 - 43

L

42 - 47

SHORT

5-7

REGULAR

7-9

**BASE OF NECK
TO TOP OF NECK**



Base of Neck



Base of Neck
to Top of Neck



Comfiwave® Kids

Optimised for infants and children with a lighter, softer feel, prioritising comfort during wear.

SHORT

REGULAR

LONG

EXTRA LONG


P

CW1-AD0S

CW1-AD0R

CW1-AD0L

CW1-AD0X



Circumferences cm

PAEDIATRIC

c WIDEST CALF 20 - 33

b ANKLE 14 - 22

Lengths cm

SHORT

REGULAR

LONG

EXTRA LONG

A-d FLOOR TO BELOW KNEE

17 - 22

22 - 27

27 - 32

32 - 37



b
Ankle



c
Widest Calf



A-d
Floor to below knee

Comfiwave[®]

Haddenham Healthcare Ltd
4D Drakes Drive
Long Crendon, HP18 9BA
United Kingdom



www.comfiwave.info



01844 208842



sales@hadhealth.com